

WARNING:

Keep out of reach of Children. In case of eye contact, rinse cautiously with water for several minutes. If eye irritation persists contact a physician.
If swallowed accidentally obtain immediate medical attention. If you obtain skin irritation contact a physician.
Flammable liquid and vapour, keep away from heat, hot surfaces, sparks, open flames and other ignition sources.

INGREDIENTS:

Active Ingredients: Quaternary ammonium compounds, benzyl-C8-18-alkyldimethyl, chlorides, Didecyldimethylammonium chloride, Amines, C12-14 alkyldimethyl N-oxides, Sodium ethylhexyl sulfate

TESTED TO EN1276, EN14476, EN13727
EN1650

BACTERIA	TEST
<i>Escherichia coli</i> (E. coli)	EN1276,
<i>Enterococcus hirae</i>	EN1276,
<i>Staphylococcus aureus</i>	EN1276,
<i>Pseudomonas aeruginosa</i>	EN1276
VIRUS	TEST
Poliiovirus	EN14476
Norovirus	EN14476
SARS-CoV-2	EN14476
Human Immunodeficiency Virus (HIV)	EN14476
Rubella Virus	EN14476
YEASTICIDAL	TEST
<i>Candida Albicans</i>	EN1650

EXPERT INFECTION CONTROL

Medisanitize

TM

**UNIVERSAL
WIPES**



ANTI-VIRAL



ANTI-BACTERIAL



ALCOHOL FREE



Effective Against:

- ✓ MRSA
- ✓ Listeria
- ✓ Salmonella
- ✓ E-Coli
- ✓ SARS-Cov-2

INSTRUCTION FOR USE:



Remove one wipe from the pack.



Working from top to bottom, wipe in an 'S' shaped pattern.



Change wipe if it becomes dry or soiled and discard. Let the surface air dry.

- ✓ EN1275
- ✓ EN1276
- ✓ EN1500
- ✓ EN1656
- ✓ EN13623
- ✓ EN13697
- ✓ EN13704
- ✓ EN13727
- ✓ EN14248
- ✓ EN14348
- ✓ EN14476
- ✓ EN1650
- ✓ ANTIBACTERIAL
- ✓ VIRUCIDAL
- ✓ 99.999% KILL RATE

MANUFACTURED BY
EXTERGEO
INDUSTRIES LTD

B5 Buckshaw Link
Buckshaw Village
Chorley
PR7 7EL
United Kingdom
Tel: 03333446019
Info@extergeo.co.uk
www.extergeo.co.uk
www.medisanitize.co.uk



PRODUCT SHEET : UNIVERSAL WIPES TUB

Medisanitize™



Description

An alcohol-free competent product contains a nonwoven fabric pre-saturated with an aqueous disinfectant solution containing a blend of high-grade biocides efficiently acting against the microbes including bacterial spores and bacteria developing resistance.

Features

- ✔ 2 in 1 use – cleaning and disinfection.
- ✔ Quick, convenient, and effective.
- ✔ Single use wipe - no cross contamination.

Packaging Details

Product	Name	Size
BT100UNI	100 Universal Tub Wipes	200mm*200mm
BT150UNI	150 Universal Tub Wipes	200mm*200mm
BT200UNI	200 Universal Tub Wipes	200mm*200mm
ST200UNI	200 Universal Probe Tub Wipes	125mm*125mm

Technical details

- ✔ ≥99.9999% kill rate against bacteria including MRSA and Mycobacteria.
- ✔ Kills at least 99.999% of Clostridium difficile spores within 5 minutes.
- ✔ 99.999% kill rate against Legionella pneumophila and Fungi.
- ✔ Yeastocidal with 99.9999% kill rate.
- ✔ Effective against all enveloped viruses including Coronavirus.

According to European standard tests:

EN 1276, EN 1656, EN 13623, EN 13697, EN 13727, EN 13704, EN 14348, EN 1650, EN 1275, EN 14476.

Applications

Excellent cleaning and disinfection of hand and multiple hard surfaces. Recommended for amateur and professional user.

Fabric Availability - Blue/White

