



BEESWIFT
FOCUSED ON SAFETY

Product Code:

CM0376

Selling UOM:

Each

Pack Size:

1

SPECIFICATION SHEET

LUXURY HOT / COLD GEL PACK LARGE RE-USABLE

BLUE DOT

Luxury large re-usable fabric Hot-Cold gel pack.

Cold therapy:

- Store in freezer for 1 hour to use for cold therapy
- Cold therapy used for swelling reduction, tooth ache, headaches and reducing fevers

Hot therapy:

- Place in microwave for 30 seconds or Boil for 7 mins to use for hot therapy
- Hot therapy used for muscle relief, back pain, sinusitis, menstrual cramps



Packaging Waste

Plastic

Carton Weight:

0.34kg

Type

Outer

Units

1

Weight(g)

2.84

Length(mm)

Width(mm)

Height(mm)